

## GABARAPL2 antibody

Cat. No. GTX31584

|              |                   |
|--------------|-------------------|
| Host         | Rabbit            |
| Clonality    | Polyclonal        |
| Isotype      | IgG               |
| Applications | WB, IHC-P, ELISA  |
| Reactivity   | Human, Mouse, Rat |

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1 - 2 µg/mL          |
| IHC-P              | 5 µg/mL              |
| ELISA              | Assay dependent      |

Not tested in other applications.

Calculated MW 14 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS                                                                                                                                                                                                                        |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| Immunogen     | GABARAPL2 antibody was raised against a 15 amino acid peptide near the carboxy terminus of human GABARAPL2. The immunogen is located within the last 50 amino acids of GABARAPL2.                                          |
| Purification  | Purified by antigen-affinity chromatography                                                                                                                                                                                |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX31584 WB Image**


WB analysis of human brain tissue lysate using GTX31584 GABARAPL2 antibody.  
Working concentration : (A) 1 and (B) 2 µg/ml


**GTX31584 IHC-P Image**

IHC-P analysis of rat brain tissue using GTX31584 GABARAPL2 antibody.  
Working concentration : 5 µg/ml


**GTX31584 IHC-P Image**

IHC-P analysis of mouse brain tissue using GTX31584 GABARAPL2 antibody.  
Working concentration : 20 µg/ml


**GTX31584 IHC-P Image**

IHC-P analysis of mouse brain tissue using GTX31584 GABARAPL2 antibody.  
Working concentration : 5 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).