

AKAP 82 blocking peptide

Cat. No. GTX31595-PEP

Application Neutralizing/Inhibition**Package**

50 µg

PRODUCT

Summary AKAP 82 peptide is used for blocking the activity of AKAP 82 antibody (catalog number GTX31595).

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution **Recommended dilution**

Neutralizing/Inhibition Assay dependent

Not tested in other applications.

PROPERTIES

Form Liquid**Buffer** PBS, 0.1% BSA**Preservative** 0.02% Sodium azide**Storage** Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.**Concentration** 200 µg/ml (Please refer to the vial label for the specific concentration.)**Conjugation** Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).