

SKI2W antibody

Cat. No. GTX31599

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 138 kDa. ([Note](#))

Product Note At least two isoforms of SKI2W are known to exist; this antibody will detect both isoforms.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	SKI2W antibody was raised against a 16 amino acid peptide near the amino terminus of human SKI2W. The immunogen is located within amino acids 60 - 110 of SKI2W.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX31599 IHC-P Image

IHC-P analysis of human stomach tissue using GTX31599 SKI2W antibody.
Working concentration : 5 µg/ml


GTX31599 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX31599 SKI2W antibody.
Dilution : 20 µg/ml


GTX31599 WB Image

WB analysis of mouse stomach tissue lysate using GTX31599 SKI2W antibody.
Working concentration : 1 µg/ml



For full product information, images and publications, please visit our [website](#).