

P2X7 blocking peptide

Cat. No. GTX31600-PEP

Applications Neutralizing /Inhibition

Package
50 µg

PRODUCT

Summary P2X7 peptide is used for blocking the activity of P2X7 antibody (catalog number GTX31600).

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

Neutralizing /Inhibition	Assay dependent
--------------------------	-----------------

Not tested in other applications.

Properties

Form Liquid

Buffer PBS, 0.1% BSA

Preservative 0.02% Sodium azide

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration 200 µg/ml (Please refer to the vial label for the specific concentration.)

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).