

KCNK1 antibody

Cat. No. GTX31601

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 38 kDa. ([Note](#))

Product Note KCNK1 antibody is predicted to not cross-react with other KCNK proteins.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KCNK1 antibody was raised against a 17 amino acid peptide near the carboxy terminus of human KCNK1. The immunogen is located within amino acids 220 - 270 of KCNK1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX31601 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX31601 KCNK1 antibody.


Working concentration : 20 µg/ml



GTX31601 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX31601 KCNK1 antibody.

Working concentration : 5 µg/ml



GTX31601 WB Image

WB analysis of 3T3 cell lysate using GTX31601 KCNK1 antibody.

Working concentration : 1 µg/ml



For full product information, images and publications, please visit our [website](#).