

ODZ3 antibody

Cat. No. GTX31616

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	301 kDa. (Note)
Product Note	This antibody is predicted to not cross-react with other members of the TENM family.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	ODZ3 / Teneurin-3 antibody was raised against a 17 amino acid peptide near the amino terminus of human ODZ3 / Teneurin-3. The immunogen is located within the first 50 amino acids of ODZ3 / Teneurin-3.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX31616 WB Image

WB analysis of 293 cell lysate using GTX31616 ODZ3 antibody.
Working concentration : 1 µg/ml



GTX31616 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX31616 ODZ3 antibody.
Working concentration : 20 µg/ml



GTX31616 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX31616 ODZ3 antibody.
Working concentration : 5 µg/ml



For full product information, images and publications, please visit our [website](#).