

Golgin 97 antibody

Cat. No. GTX31624

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 88 kDa. ([Note](#))

Product Note At least two isoforms of GOLGIN97 are known to exist; this antibody will detect only the longer isoform. GOLGIN97 antibody is predicted to not cross-react with other GOLGIN family proteins.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	GOLGIN97 antibody was raised against a 17 amino acid peptide near the amino terminus of human GOLGIN97. The immunogen is located within amino acids 30 - 80 of GOLGIN97.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX31624 IHC-P Image


IHC-P analysis of rat stomach tissue using GTX31624 Golgin 97 antibody.

Working concentration : 20 µg/ml


GTX31624 IHC-P Image

IHC-P analysis of rat stomach tissue using GTX31624 Golgin 97 antibody.

Working concentration : 5 µg/ml


GTX31624 WB Image

WB analysis of human stomach tissue lysate in (A) the absence and (B) the presence of blocking peptide using GTX31624 Golgin 97 antibody.

Working concentration : 1 µg/ml



For full product information, images and publications, please visit our [website](#).