

Livin antibody

Cat. No. GTX31663

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	33 kDa. (Note)
Product Note	A lower but much weaker band at 30 kDa was detected in Raji cell lysate, which may represent the short form of Livin.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Livin antibody was raised against a 17 amino acid peptide near the carboxy terminus of human Livin. The immunogen is located within the last 50 amino acids of Livin.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx31663 WB Image


WB analysis of Raji cell lysate using GTx31663 Livin antibody.

Working concentration : 0.5 µg/ml


GTx31663 IHC-P Image

IHC-P analysis of human small intestine tissue using GTx31663 Livin antibody.

Working concentration : 5 µg/ml


GTx31663 IHC-P Image

IHC-P analysis of human small intestine tissue using GTx31663 Livin antibody.

Working concentration : 20 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).