

TIRAP antibody

Cat. No. GTX31676

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	2 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 26 kDa. ([Note](#))

Product Note At least three isoforms of TIRAP are known to exist; this antibody detects all three isoforms.

Properties


Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	TIRAP antibody was raised against a peptide corresponding to 15 amino acids near the C-terminus of mouse TIRAP. The immunogen is located within the last 50 amino acids of TIRAP.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX31676 IHC-P Image

IHC-P analysis of human heart tissue using GTX31676 TIRAP antibody.
Working concentration : 2 µg/ml


GTX31676 IHC-P Image

IHC-P analysis of human heart tissue using GTX31676 TIRAP antibody.
Working concentration : 10 µg/ml


GTX31676 WB Image

WB analysis of human heart tissue lysate using GTX31676 TIRAP antibody.
Working concentration : (A) 1 and (B) 2 µg/ml



For full product information, images and publications, please visit our [website](#).