

TLR11 antibody

Cat. No. GTX31682

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, ELISA
Reactivity	Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5 - 1 µg/mL
ICC/IF	2 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 106 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	TLR11 antibody was raised against a peptide corresponding to 15 amino acids near the carboxy terminus of mouse TLR11. The immunogen is located within the last 50 amino acids of TLR11.
Purification	Purified by ion exchange chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

DATA IMAGES

GTX31682 ICC/IF Image

ICC/IF analysis of RAW264.7 cells using GTX31682 TLR11 antibody.
Working concentration : 2 µg/ml


GTX31682 ICC/IF Image

ICC/IF analysis of RAW264.7 cells using GTX31682 TLR11 antibody.
Working concentration : 10 µg/ml


GTX31682 WB Image

WB analysis of RAW264.7 cell lysates using GTX31682 TLR11 antibody.
Dilution : 0.5µg/ml (Lane A) and 1 µg/ml (Lane B)



For full product information, images and publications, please visit our [website](https://www.genetex.com).