

# LSD1 blocking peptide

**Cat. No. GTX31730-PEP**

## Applications

Neutralizing /Inhibition

## Package

50 µg

## PRODUCT

### Summary

LSD1 peptide is used for blocking the activity of LSD1 antibody (catalog number GTX31730).

### Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

### Suggested dilution

### Recommended dilution

Neutralizing /Inhibition

Assay dependent

Not tested in other applications.

## Properties

**Form** Liquid

**Buffer** PBS, 0.1% BSA

**Preservative** 0.02% Sodium azide

**Storage** Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

**Concentration** 200 µg/ml (Please refer to the vial label for the specific concentration.)

**Conjugation** Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).