

## Bag1 antibody

Cat. No. GTX31737

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Application	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

**Package**  
100 µg

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	2 µg/mL
ELISA	Assay dependent

Not tested in other applications.

**Calculated MW** 39 kDa. ( [Note](#) )

## PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	BAG-1 antibody was raised against a 14 amino acid synthetic peptide from near the amino terminus of human BAG-1. The immunogen is located within the first 50 amino acids of BAG-1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX31737 IHC-P Image**

IHC-P analysis of human lymph node tissue using GTX31737 Bag1 antibody.  
Working concentration : 2 µg/ml


**GTX31737 WB Image**

WB analysis of PC-3 cell lysate using GTX31737 Bag1 antibody.  
Working concentration : (A) 1 and (B) 2 µg/ml


**GTX31737 IHC-P Image**

IHC-P analysis of human lymph node tissue using GTX31737 Bag1 antibody.  
Working concentration : 20 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).