

AS160 / TBC1D4 blocking peptide

Cat. No. GTX31767-PEP

Applications

Neutralizing /Inhibition

Package

50 µg

PRODUCT

Summary

AS160 / TBC1D4 peptide is used for blocking the activity of AS160 / TBC1D4 antibody (catalog number GTX31767). It usually blocks the antibody activity completely in Western blot by incubating the peptide with equal volume of antibody for 30 min at 37°C.

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

Recommended dilution

Neutralizing /Inhibition

Assay dependent

Not tested in other applications.

Properties

Form

Liquid

Buffer

PBS, 0.1% BSA

Preservative

0.02% Sodium azide

Storage

Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration

200 µg/ml (Please refer to the vial label for the specific concentration.)

Conjugation

Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).