

Alpha fetoprotein / AFP blocking peptide

Cat. No. GTX31809-PEP

Applications

Neutralizing /Inhibition

Package

50 µg

PRODUCT

Summary

Alpha fetoprotein / AFP peptide is used for blocking the activity of Alpha fetoprotein / AFP antibody (catalog number GTX31809).

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

Recommended dilution

Neutralizing /Inhibition

Assay dependent

Not tested in other applications.

Properties

Form

Liquid

Buffer

PBS, 0.1% BSA

Preservative

0.02% Sodium azide

Storage

Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration

200 µg/ml (Please refer to the vial label for the specific concentration.)

Conjugation

Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx31809-PEP Image

WB analysis of fetal human liver tissue lysate using GTx31809 Alpha fetoprotein / AFP antibody. The lane on the right was treated with GTx31809-PEP Alpha fetoprotein / AFP blocking peptide.

Loading : 10 μ g

Dilution : 0.5 μ g/mL



For full product information, images and publications, please visit our [website](#).