

SLC39A9 antibody

Cat. No. GTX31817

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse

References (3)

★★★★★ Review (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	2.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 32 kDa. ([Note](#))

Product Note ZIP9 antibody is predicted to not cross-react with other ZIP family members.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	SLC39A9 antibody was raised against a 18 amino acid synthetic peptide near the amino terminus of human SLC39A9. The immunogen is located within amino acids 70 - 120 of SLC39A9.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX31817 IHC-P Image

IHC-P analysis of human brain tissue using GTX31817 SLC39A9 antibody.
Working concentration : 2.5 µg/ml



GTX31817 WB Image

WB analysis of HepG2 cell lysate using GTX31817 SLC39A9 antibody.
Working concentration : (A) 1 and (B) 2 µg/ml



GTX31817 IHC-P Image

IHC-P analysis of human brain tissue using GTX31817 SLC39A9 antibody.
Working concentration : 20 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).