

NELF antibody

Cat. No. GTX31872

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 60 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rabbit polyclonal NELF antibody was raised against a 17 amino acid peptide near the center of human NELF. The immunogen is located within amino acids 230 - 280 of NELF.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX31872 IHC-P Image

IHC-P analysis of rat brain tissue using GTX31872 NELF antibody.
Working concentration : 20 µg/ml




GTX31872 WB Image

WB analysis of HeLa cell lysate using GTX31872 NELF antibody.
Working concentration : (A) 1 and (B) 2 µg/ml



GTX31872 IHC-P Image

IHC-P analysis of rat brain tissue using GTX31872 NELF antibody.
Working concentration : 2.5 µg/ml



GTX31872 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX31872 NELF antibody.
Working concentration : 5 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).