

KREMEN1 antibody

Cat. No. GTX31876

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.125 - 0.25 µg/mL
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 52 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rabbit polyclonal KREMEN1 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human KREMEN1. The immunogen is located within the last 50 amino acids of KREMEN1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX31876 IHC-P Image

IHC-P analysis of human small intestine tissue using GTX31876 KREMEN1 antibody.
Working concentration : 20 µg/ml



GTX31876 IHC-P Image

IHC-P analysis of human small intestine tissue using GTX31876 KREMEN1 antibody.
Working concentration : 5 µg/ml



GTX31876 WB Image

WB analysis of rat small intestine tissue lysate using GTX31876 KREMEN1 antibody.
Working concentration : (A) 0.125 and (B) 0.25 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).