

TFPI antibody

Cat. No. GTX31885

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	35 kDa. (Note)
Product Note	At least two isoforms of TFPI are known to exist; this antibody will detect both isoforms.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	TFPI antibody was raised against a 13 amino acid peptide near the amino terminus of human TFPI. The immunogen is located within amino acids 70 - 120 of TFPI.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 09 / 14 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX31885 IHC-P Image

IHC-P analysis of rat small intestine tissue using GTX31885 TFPI antibody.


Working concentration : 20 µg/ml



GTX31885 IHC-P Image

IHC-P analysis of rat small intestine tissue using GTX31885 TFPI antibody.

Working concentration : 5 µg/ml



GTX31885 WB Image

WB analysis of rat small intestine tissue lysate using GTX31885 TFPI antibody.

Working concentration : 1 µg/ml



For full product information, images and publications, please visit our [website](#).