

SPNS1 antibody

Cat. No. GTX31907

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	2.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 57 kDa. ([Note](#))

Product Note This antibody is human specific. At least four isoforms are known to exist. This antibody is predicted to not cross-react with other members of the spinster family of proteins.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	SPNS1 antibody was raised against an 18 amino acid peptide near the amino terminus of human SPNS1. The immunogen is located within the first 50 amino acids of SPNS1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX31907 WB Image

WB analysis of HepG2 cell lysate using GTX31907 SPNS1 antibody.
Working concentration : (A) 1 and (B) 2 µg/ml



GTX31907 IHC-P Image

IHC-P analysis of human lung tissue using GTX31907 SPNS1 antibody.
Working concentration : 20 µg/ml



GTX31907 IHC-P Image

IHC-P analysis of human lung tissue using GTX31907 SPNS1 antibody.
Working concentration : 2.5 µg/ml



For full product information, images and publications, please visit our [website](#).