

CADM3 antibody

Cat. No. GTX31935

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 43 kDa. ([Note](#))

Product Note At least three isoforms are known to exist; this antibody will only detect the two longest isoforms. CADM3 antibody is predicted to not cross-react with other members of the CADM protein family.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	CADM3 antibody was raised against a 16 amino acid peptide near the center of human CADM3. The immunogen is located within amino acids 220 - 270 of CADM3.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX31935 WB Image

WB analysis of human brain tissue lysate using GTX31935 CADM3 antibody.

Dilution : 1 µg/ml (Lane A) and 2 µg/ml (Lane B)



GTX31935 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX31935 CADM3 antibody.

Working concentration : 20 µg/ml



GTX31935 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX31935 CADM3 antibody.

Working concentration : 5 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).