

JAK2 blocking peptide

Cat. No. GTX31943-PEP

Applications Neutralizing /Inhibition

Package
50 µg

PRODUCT

Summary JAK2 peptide is used for blocking the activity of JAK2 antibody (catalog number GTX31943).

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Neutralizing /Inhibition	Assay dependent
Not tested in other applications.	

Properties

Form Liquid

Buffer PBS, 0.1% BSA

Preservative 0.02% Sodium azide

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration 200 µg/ml (Please refer to the vial label for the specific concentration.)

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).