

Adenosine A2b Receptor blocking peptide

Cat. No. GTX31945-PEP

Applications

Neutralizing /Inhibition

Package

50 µg

PRODUCT

Summary

Adenosine A2b Receptor peptide is used for blocking the activity of Adenosine A2b Receptor antibody (catalog number GTX31945).

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

Recommended dilution

Neutralizing /Inhibition

Assay dependent

Not tested in other applications.

Properties

Form Liquid

Buffer PBS, 0.1% BSA

Preservative 0.02% Sodium azide

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration 200 µg/ml (Please refer to the vial label for the specific concentration.)

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).