

SFRS3 antibody

Cat. No. GTX31947

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
ICC/IF	2 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 19 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	SFRS3 antibody was raised against a 19 amino acid peptide near the center of human SFRS3. The immunogen is located within amino acids 80 - 130 of SFRS3.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX31947 WB Image

WB analysis of K562 cell lysate using GTX31947 SFRS3 antibody.
Working concentration : (A) 1 and (B) 2 µg/ml



GTX31947 ICC/IF Image

ICC/IF analysis of K562 cells using GTX31947 SFRS3 antibody.
Working concentration : 2 µg/ml



GTX31947 ICC/IF Image

ICC/IF analysis of K562 cells using GTX31947 SFRS3 antibody.
Working concentration : 2 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).