

RFX1 antibody

Cat. No. GTX31948

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
ICC/IF	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 105 kDa. ([Note](#))

Product Note RFX1 antibody is human specific. At least two isoforms of RFX1 are known to exist; this antibody will detect both isoforms.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	RFX1 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human RFX1. The immunogen is located within the last 50 amino acids of RFX1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 07 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTx31948 WB Image

WB analysis of HeLa cell lysate using GTx31948 RFX1 antibody.


Working concentration : (A) 1 and (B) 2 µg/ml



GTx31948 ICC/IF Image

ICC/IF analysis of HeLa cells using GTx31948 RFX1 antibody.

Working concentration : 5 µg/ml



GTx31948 ICC/IF Image

ICC/IF analysis of HeLa cells using GTx31948 RFX1 antibody.

Green : Primary antibody

Blue : DAPI

Dilution : 20 µg/ml



For full product information, images and publications, please visit our [website](#).