

ZMYM1 antibody

Cat. No. GTX32000

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 129 kDa. ([Note](#))

Product Note At least two isoforms of ZMYM1 are known to exist; this antibody will detect only the shorter isoform. ZMYM1 antibody is predicted to not cross-react with other ZMYM protein family members.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	ZMYM1 antibody was raised against a 17 amino acid synthetic peptide near the carboxy terminus of human ZMYM1 .The immunogen is located within amino acids 990 - 1040 of ZMYM1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 07 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX32000 WB Image

WB analysis of rat lung tissue lysate in (A) the absence and (B) the presence of blocking peptide using GTX32000 ZMYM1 antibody.

Working concentration : 1 µg/ml



GTX32000 IHC-P Image

IHC-P analysis of rat lung tissue using GTX32000 ZMYM1 antibody.

Working concentration : 5 µg/ml



For full product information, images and publications, please visit our [website](#).