

PION antibody

Cat. No. GTX32001

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.25 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	98 kDa. (Note)
Product Note	Multiple isoforms of PION are known to exist. PION antibody is predicted to not cross-react with other F-box protein family members.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	PION antibody was raised against a 19 amino acid synthetic peptide near the carboxy terminus of human PION. The immunogen is located within amino acids 770 - 820 of PION.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTx32001 IHC-P Image

IHC-P analysis of human brain tissue using GTx32001 PION antibody.
Working concentration : 5 µg/ml



GTx32001 IHC-P Image

IHC-P analysis of human brain tissue using GTx32001 PION antibody.
Working concentration : 20 µg/ml



GTx32001 IHC-P Image

IHC-P analysis of mouse brain tissue using GTx32001 PION antibody.
Working concentration : 5 µg/ml



For full product information, images and publications, please visit our [website](#).