

VPS39 antibody

Cat. No. GTX32005

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 102 kDa. ([Note](#))

Product Note At least four isoforms of VPS39 are known to exist; this antibody will detect all four. VPS39 antibody is predicted not to cross-react with TRAP1.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	VPS39 antibody was raised against a 18 amino acid synthetic peptide near the carboxy terminus of human VPS39. The immunogen is located within amino acids 710 - 760 of VPS39.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTx32005 WB Image

WB analysis of rat liver tissue lysate in (A) the absence and (B) the presence of blocking peptide using GTx32005 VPS39 antibody.

Working concentration : 0.5 µg/ml



GTx32005 IHC-P Image


IHC-P analysis of human liver tissue using GTx32005 VPS39 antibody.

Working concentration : 5 µg/ml

Green : Primary antibody

Blue : DAPI

Red : Actin



GTx32005 IHC-P Image

IHC-P analysis of mouse liver tissue using GTx32005 VPS39 antibody.

Working concentration : 20 µg/ml



For full product information, images and publications, please visit our [website](#).