

KAL1 antibody

Cat. No. GTX32013

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
ICC/IF	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 76 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KAL1 antibody was raised against an 18 amino acid synthetic peptide near the center of human KAL1. The immunogen is located within amino acids 430 - 480 of KAL1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

DATA IMAGES




GTX32013 ICC/IF Image

ICC/IF analysis of MCF7 cells using GTX32013 KAL1 antibody.
Working concentration : 20 µg/ml



GTX32013 WB Image

WB analysis of MCF7 cell lysate using GTX32013 KAL1 antibody.
Working concentration : 1 µg/ml



GTX32013 ICC/IF Image

ICC/IF analysis of MCF7 cells using GTX32013 KAL1 antibody.
Working concentration : 5 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).