

HDGFRP2 antibody

Cat. No. GTX32024

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5 - 1 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 74 kDa. ([Note](#))

Product Note At least two isoforms of HDGFRP2 are known to exist; this antibody will recognize both isoforms. This antibody is predicted to not cross-react with other HDGFRP family members

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	HDGFRP2 antibody was raised against an 18 amino acid synthetic peptide near the amino terminus of human HDGFRP2. The immunogen is located within amino acids 50 - 100 of HDGFRP2.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX32024 IHC-P Image

IHC-P analysis of rat liver tissue using GTX32024 HDGFRP2 antibody.


Working concentration : 20 µg/ml



GTX32024 IHC-P Image

IHC-P analysis of rat liver tissue using GTX32024 HDGFRP2 antibody.

Working concentration : 5 µg/ml



GTX32024 WB Image

WB analysis of rat lung tissue lysate using GTX32024 HDGFRP2 antibody.

Working concentration : 2 µg/ml



For full product information, images and publications, please visit our [website](#).