

SPINK2 antibody

Cat. No. GTX32051

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse

References (2)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
IHC-P	2.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 10 kDa. ([Note](#))

Product Note SPINK2 antibody is predicted to not cross-react with other SPINK family members.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	SPINK2 antibody was raised against a 16 amino acid synthetic peptide near the carboxy terminus of mouse SPINK2. The immunogen is located within the last 50 amino acids of SPINK2.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX32051 IHC-P Image

IHC-P analysis of mouse testis tissue using GTX32051 SPINK2 antibody.


Working concentration : 20 µg/ml



GTX32051 IHC-P Image

IHC-P analysis of mouse testis tissue using GTX32051 SPINK2 antibody.

Working concentration : 2.5 µg/ml



GTX32051 WB Image

WB analysis of human testis tissue lysate in (A) the absence and (B) the presence of blocking peptide using GTX32051 SPINK2 antibody.

Working concentration : 1 µg/ml



For full product information, images and publications, please visit our [website](#).