

MLIP antibody

Cat. No. GTX32056

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	50 kDa. (Note)
Product Note	At least three isoforms of MLIP are known to exist.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	MLIP antibody was raised against a 16 amino acid synthetic peptide near the carboxy terminus of human MLIP. The immunogen is located within the last 50 amino acids of MLIP.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTx32056 WB Image

WB analysis of 293 cell lysate in (A) the absence and (B) the presence of blocking peptide using GTx32056 MLIP antibody.


Working concentration : 1 µg/ml



GTx32056 IHC-P Image

IHC-P analysis of human heart tissue using GTx32056 MLIP antibody.

Working concentration : 5 µg/ml



GTx32056 IHC-P Image

IHC-P analysis of human heart tissue using GTx32056 MLIP antibody.

Working concentration : 20 µg/ml



For full product information, images and publications, please visit our [website](#).