

GZF1 antibody

Cat. No. GTX32070

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
IHC-P	2.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 80 kDa. ([Note](#))

Product Note GZF1 antibody is human specific. At least four isoforms of GZF1 are known to exist.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	GZF1 antibody was raised against an 18 amino acid synthetic peptide near the carboxy terminus of human GZF1. The immunogen is located within amino acids 600 - 650 of GZF1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX32070 WB Image

WB analysis of human heart tissue lysate using GTX32070 GZF1 antibody.


Working concentration : 1 µg/ml



GTX32070 IHC-P Image

IHC-P analysis of human kidney tissue using GTX32070 GZF1 antibody.

Working concentration : 20 µg/ml



GTX32070 IHC-P Image

IHC-P analysis of human kidney tissue using GTX32070 GZF1 antibody.

Working concentration : 2.5 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).