

LYZL1 blocking peptide

Cat. No. GTX32083-PEP**Applications**

Neutralizing /Inhibition

Package

50 µg

PRODUCT**Summary**

LYZL1 peptide is used for blocking the activity of LYZL1 antibody (catalog number GTX32083).

Applications**Application Note**

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution**Recommended dilution**

Neutralizing /Inhibition

Assay dependent

Not tested in other applications.

Properties**Form**

Liquid

Buffer

PBS, 0.1% BSA

Preservative

0.02% Sodium azide

Storage

Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration

200 µg/ml (Please refer to the vial label for the specific concentration.)

Conjugation

Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).