

MOSPD1 antibody

Cat. No. GTX32111

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 24 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	MOSPD1 antibody was raised against a 19 amino acid peptide near the center of human MOSPD1 .The immunogen is located within amino acids 110 - 160 of MOSPD1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX32111 IHC-P Image

IHC-P analysis of mouse testis tissue using GTX32111 MOSPD1 antibody.
Working concentration : 2.5 µg/ml


GTX32111 IHC-P Image

IHC-P analysis of mouse testis tissue using GTX32111 MOSPD1 antibody.
Working concentration : 20 µg/ml


GTX32111 WB Image

WB analysis of human brain tissue lysate in (A) the absence and (B) the presence of blocking peptide using GTX32111 MOSPD1 antibody.
Working concentration : 1 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).