

MTERFD3 antibody

Cat. No. GTX32114

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	44 kDa. (Note)
Product Note	MTERFD3 antibody is predicted to not cross-react with other MTERFD protein family members.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	MTERFD3 antibody was raised against an 18 amino acid peptide near the amino terminus of human MTERFD3. The immunogen is located within amino acids 20 - 70 of MTERFD3.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX32114 WB Image

WB analysis of human testis tissue lysate using GTX32114 MTERFD3 antibody.

Working concentration : 1 µg/ml



GTX32114 IHC-P Image

IHC-P analysis of mouse testis tissue using GTX32114 MTERFD3 antibody.

Working concentration : 2.5 µg/ml



For full product information, images and publications, please visit our [website](#).