

LSD2 / AOF1 antibody

Cat. No. GTX32122

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Human, Mouse

References (2)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 92 kDa. ([Note](#))

Product Note At least two isoforms of KDM1B are known to exist; this KDM1B antibody will detect both isoforms.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KDM1B antibody was raised against a 17 amino acid peptide near the amino terminus of human KDM1B. The immunogen is located within amino acids 80 - 130 of KDM1B.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX32122 WB Image

WB analysis of K562 cell lysate in (A) the absence and (B) the presence of blocking peptide using GTX32122 KDM1B / AOF1 antibody.

Working concentration : 0.5 µg/ml



For full product information, images and publications, please visit our [website](#).