

TRIP12 antibody

Cat. No. GTX32138

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
ICC/IF	Assay dependent
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 220 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	TRIP12 antibody was raised against a 18 amino acid peptide near the carboxy terminus of human TRIP12. The immunogen is located within amino acids 1930 - 1980 of TRIP12.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX32138 IHC-P Image


IHC-P analysis of rat colon tissue using GTX32138 TRIP12 antibody.
Working concentration : 5 µg/ml


GTX32138 IHC-P Image

IHC-P analysis of rat colon tissue using GTX32138 TRIP12 antibody.
Working concentration : 20 µg/ml


GTX32138 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX32138 TRIP12 antibody.
Dilution : 2 µg/ml


GTX32138 WB Image

WB analysis of rat colon tissue lysate using GTX32138 TRIP12 antibody.
Working concentration : 1 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).