

ULBP4 antibody

Cat. No. GTX32159

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, ELISA
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
ICC/IF	2.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	30 kDa. (Note)
Product Note	At least four isoforms are known to exist; this antibody will detect all four.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	ULBP4 antibody was raised against an 18 amino acid peptide near the center of human ULBP4. The immunogen is located within amino acids 160 - 210 of ULBP4.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTx32159 WB Image

WB analysis of EL4 cell lysate using GTx32159 ULBP4 antibody.

Working concentration : 1 µg/ml



GTx32159 ICC/IF Image

ICC/IF analysis of K562 cells using GTx32159 ULBP4 antibody.

Working concentration : 2.5 µg/ml



For full product information, images and publications, please visit our [website](#).