

# PLCB3 (phospho Ser537) antibody

**Cat. No. GTX32237**

|                     |                   |
|---------------------|-------------------|
| <b>Host</b>         | Rabbit            |
| <b>Clonality</b>    | Polyclonal        |
| <b>Isotype</b>      | IgG               |
| <b>Applications</b> | WB                |
| <b>Reactivity</b>   | Human, Mouse, Rat |

**Package**  
100 µl

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500 - 1:1000       |

Not tested in other applications.

**Calculated MW** 139 kDa. ([Note](#))

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | 0.42% Potassium Phosphate, 0.87% NaCl, 30% Glycerol                                                                                                                                                                        |
| <b>Preservative</b>  | 0.01% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| <b>Immunogen</b>     | KLH-conjugated synthetic peptide encompassing a sequence within the center region of PLCB3. The exact sequence is proprietary.                                                                                             |
| <b>Purification</b>  | Purified by antigen-affinity chromatography                                                                                                                                                                                |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX32237 WB Image

WB analysis of PMA-treated HeLa (A), UV-treated MCF7 (B), mouse brain (C), rat brain (D) whole cell lysates using GTX32237 PLCB3 (phospho Ser537) antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).