

## NFkB p105 (phospho Ser932) antibody

Cat. No. GTX32341

|              |                   |
|--------------|-------------------|
| Host         | Rabbit            |
| Clonality    | Polyclonal        |
| Isotype      | IgG               |
| Applications | WB, IHC-P         |
| Reactivity   | Human, Mouse, Rat |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500 - 1:1000       |
| IHC-P              | 1:100 - 1:200        |

Not tested in other applications.

Calculated MW 105 kDa. ([Note](#))

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | 0.42% Potassium Phosphate, 0.87% NaCl, 30% Glycerol                                                                                                                                                                        |
| Preservative  | 0.01% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of NFkB p105. The exact sequence is proprietary.                                                                                         |
| Purification  | Purified by antigen-affinity chromatography                                                                                                                                                                                |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTx32341 WB Image

WB analysis of LPS-treated HeLa (A), LPS-treated NIH3T3 (B), LPS-treated Raw264.7 (C), rat kidney (D) whole cell lysates using GTx32341 NFκB p105 (phospho Ser932) antibody.



### GTx32341 IHC-P Image

IHC-P analysis of formalin fixed human breast cancer tissue section using GTx32341 NFκB p105 (phospho Ser932) antibody.

Antigen retrieval : Heat mediated antigen retrieval with sodium citrate buffer (pH 6.0)



For full product information, images and publications, please visit our [website](https://www.genetex.com).