

Cytokeratin 17 antibody

Cat. No. GTX32548

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:20 - 1:100

Not tested in other applications.

Calculated MW 48 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 133-432 of human KRT17 (NP_000413.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX32548 WB Image

WB analysis of various sample lysates using GTX32548 Cytokeratin 17 antibody.


Dilution : 1:1000

Loading : 25µg per lane


GTX32548 ICC/IF Image

ICC/IF analysis of A549 cells using GTX32548 Cytokeratin 17 antibody.

Blue : DAPI


GTX32548 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX32548 Cytokeratin 17 antibody.

Blue : DAPI

Dilution : 1:200



For full product information, images and publications, please visit our [website](https://www.genetex.com).