

HGF antibody

Cat. No. GTX32651

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	1:50 - 1:100

Not tested in other applications.

Calculated MW 83 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-285 of human HGF (NP_001010933.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx32651 WB Image

WB analysis of HepG2 cell lysate using GTx32651 HGF antibody.


Dilution : 1:500

Loading : 25µg per lane


GTx32651 IHC-P Image

IHC-P analysis of rat kidney tissue using GTx32651 HGF antibody.

Dilution : 1:100


GTx32651 IHC-P Image

IHC-P analysis of mouse spinal cord tissue using GTx32651 HGF antibody.

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).