

# K-Ras antibody

**Cat. No. GTX32694**

<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	WB, ICC/IF
<b>Reactivity</b>	Human, Mouse, Rat

**Package**  
100 µl

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:1000
ICC/IF	1:50 - 1:200

Not tested in other applications.

**Calculated MW** 22 kDa. ( [Note](#) )

**Product Note** This antibody may cross react with N-Ras and H-Ras proteins.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 50% Glycerol
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Recombinant protein of human KRAS
<b>Purification</b>	Purified by affinity chromatography
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


**DATA IMAGES**

**GTX32694 ICC/IF Image**

ICC/IF analysis of HeLa cells using GTX32694 K-Ras antibody.

Blue : DAPI


Dilution : 1:100


**GTX32694 WB Image**

WB analysis of various sample lysates using GTX32694 K-Ras antibody.

Dilution : 1:1000

Loading : 25µg per lane


**GTX32694 ICC/IF Image**

ICC/IF analysis of NIH/3T3 cells using GTX32694 K-Ras antibody.

Blue : DAPI

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).