

Prohibitin 2 antibody

Cat. No. GTX32812

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

References (2)

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:50 - 1:100
IHC-P	1:50 - 1:100

Not tested in other applications.

Calculated MW 33 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-299 of human PHB2 (NP_001138303.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX32812 IHC-P Image

IHC-P analysis of human esophagus tissue using GTX32812 Prohibitin 2 antibody.
Dilution : 1:100


GTX32812 ICC/IF Image

ICC/IF analysis of L929 cells using GTX32812 Prohibitin 2 antibody.
Blue : DAPI
Dilution : 1:100


GTX32812 ICC/IF Image

ICC/IF analysis of U2OS cells using GTX32812 Prohibitin 2 antibody.
Blue : DAPI
Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).