

delta Sarcoglycan antibody

Cat. No. GTX32871

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

References (2)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
IHC-P	1:50 - 1:100

Not tested in other applications.

Calculated MW	32 kDa. (Note)
Product Note	KO/KD validation is based on published data (PMID: 31689918).

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 61-290 of human SGCD (NP_000328.2).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX32871 WB Image

WB analysis of various sample lysates using GTX32871 delta Sarcoglycan antibody.

Dilution : 1:1000


Loading : 25µg per lane



GTX32871 IHC-P Image

IHC-P analysis of rat heart tissue using GTX32871 delta Sarcoglycan antibody.

Dilution : 1:100



GTX32871 IHC-P Image

The data was published in the journal Int J Mol Sci in 2019. [PMID: 31689918](https://pubmed.ncbi.nlm.nih.gov/31689918/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).