

Anterior Gradient 2 antibody

Cat. No. GTX32995

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:50 - 1:200
IHC-P	1:50 - 1:200

Not tested in other applications.

Calculated MW 20 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-175 of human AGR2 (NP_006399.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX32995 IHC-P Image

IHC-P analysis of mouse lung tissue using GTX32995 Anterior Gradient 2 antibody.
Dilution : 1:100



GTX32995 ICC/IF Image

ICC/IF analysis of L929 cells using GTX32995 Anterior Gradient 2 antibody.
Blue : DAPI
Dilution : 1:100



GTX32995 ICC/IF Image

ICC/IF analysis of U2OS cells using GTX32995 Anterior Gradient 2 antibody.
Blue : DAPI
Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).