

Lamin B2 antibody

Cat. No. GTX33292

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Application	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:50 - 1:100
IHC-P	1:50 - 1:200

Not tested in other applications.

Calculated MW 70 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 341-620 of human LMNB2 (NP_116126.3).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES


GTX33292 IHC-P Image


IHC-P analysis of mouse brain tissue using GTX33292 Lamin B2 antibody.

Dilution : 1:100


GTX33292 IHC-P Image


IHC-P analysis of human adenomyosis tissue using GTX33292 Lamin B2 antibody.

Dilution : 1:100


GTX33292 ICC/IF Image

ICC/IF analysis of U2OS cells using GTX33292 Lamin B2 antibody.

Blue : DAPI


GTX33292 IHC-P Image

IHC-P analysis of rat brain tissue using GTX33292 Lamin B2 antibody.

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).