

## OAT antibody

Cat. No. GTX33373

|                     |                   |
|---------------------|-------------------|
| <b>Host</b>         | Rabbit            |
| <b>Clonality</b>    | Polyclonal        |
| <b>Isotype</b>      | IgG               |
| <b>Applications</b> | WB, IHC-P, IP     |
| <b>Reactivity</b>   | Human, Mouse, Rat |

**Package**  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500 - 1:2000       |
| IHC-P              | 1:50 - 1:200         |
| IP                 | 1:50 - 1:100         |

Not tested in other applications.

**Calculated MW** 49 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS, 50% Glycerol                                                                                                                                                                                                          |
| <b>Preservative</b>  | 0.02% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| <b>Immunogen</b>     | Recombinant fusion protein containing a sequence corresponding to amino acids 160-439 of human OAT (NP_000265.1).                                                                                                          |
| <b>Purification</b>  | Purified by affinity chromatography                                                                                                                                                                                        |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTX33373 IP Image

IP analysis of MCF-7 cell lysate using GTX33373 OAT antibody.

Antibody amount : 3 $\mu$ g / 200 $\mu$ g lysate


Dilution : 1:1000



### GTX33373 IHC-P Image

IHC-P analysis of human esophageal cancer tissue using GTX33373 OAT antibody.


Dilution : 1:150



### GTX33373 IHC-P Image

IHC-P analysis of rat kidney tissue using GTX33373 OAT antibody.

Dilution : 1:150



### GTX33373 WB Image

WB analysis of various sample lysates using GTX33373 OAT antibody.

Dilution : 1:1000

Loading : 25 $\mu$ g per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).