

PPP4R1 antibody

Cat. No. GTX33426

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
IHC-P	1:50 - 1:100
IP	1:50 - 1:100

Not tested in other applications.

Calculated MW 107 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 324-647 of human PPP4R1 (NP_001035847.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX33426 IHC-P Image


IHC-P analysis of human liver tissue using GTX33426 PPP4R1 antibody.
Dilution : 1:100


GTX33426 WB Image

WB analysis of various sample lysates using GTX33426 PPP4R1 antibody.
Dilution : 1:1000
Loading : 25µg per lane


GTX33426 IHC-P Image

IHC-P analysis of human stomach tissue using GTX33426 PPP4R1 antibody.
Dilution : 1:100


GTX33426 IP Image

IP analysis of 293T cell lysate using GTX33426 PPP4R1 antibody.
Antibody amount : 3µg / 200µg lysate
Dilution : 1:1000



For full product information, images and publications, please visit our [website](http://www.genetex.com).