

RPS10 antibody

Cat. No. GTX33475

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IP
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:20 - 1:50
IHC-P	1:50 - 1:200
IP	1:20 - 1:50

Not tested in other applications.

Calculated MW 19 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.05% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human RPS10 (NP_001190174.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX33475 IP Image

IP analysis of MCF-7 cell lysate using GTX33475 RPS10 antibody.

Antibody amount : 1 μ g / 200 μ g lysate


Dilution : 1:1000


GTX33475 WB Image

WB analysis of various sample lysates using GTX33475 RPS10 antibody.

Dilution : 1:1000

Loading : 25 μ g per lane


GTX33475 ICC/IF Image

ICC/IF analysis of MCF7 cells using GTX33475 RPS10 antibody.

Blue : DAPI



For full product information, images and publications, please visit our [website](https://www.genetex.com).